


REFRIGERATOR TRAY ASSEMBLY INSTRUCTIONS

This tray was designed to be attached to wire grid type refrigerator shelves.

Tools required: a 7/16" wrench or socket and a wide blade screwdriver.

Using the enclosed sketch as reference, hold the tray under a shelf and as close to the front edge of the shelf as possible, without interfering with the tray's lock.

Lay one of the two brackets across the top of the shelf so that the bracket holes line up with those in the top of the tray and then drop two of the enclosed bolts through each pair of holes and tighten nuts with fingers on each of these bolts from the bottom of the tray. Repeat this process with the second bracket. If possible, this process is much easier if you can remove the refrigerator shelf unit on which the tray is to be installed.



If you are mounting to a glass shelf or location other than a refrigerator shelf, the double sided tape included with your box can be used.

Put the wrench or socket on each nut from the bottom of the tray and tighten with the screwdriver from the top. Do not completely tighten all the bolts/nuts until you have determined the final position. Note: these are special nuts that require a fair amount of effort to install...however, once installed they can not be loosened without access to the bottom of the tray.

Turn the key in the lock so that the lock tab moves out of the way. Slide in the plastic box and turn the key so that the key lock tab is pointing down and is in front of the box (this is a visual indication that the box is secured within the locking tray).

If you have any questions about the installation and/or use of this locking tray, please call our toll free number for assistance (1-800-848-1633).



Call Free 1 800 848 1633 or 1 888 HCL-INTL • Fax Free 1 800 447 2923
Web Site: www.HealthCareLogistics.com • E-mail: hcl@HealthCareLogistics.com