

Instructions for the Arm and Leg Stand (#10656)



1. The Arm and Leg Stand should include:

- A •1- tripod
- B •1-5-footed plastic mid-section
- C •1-clamp disk
- D •2-feet with locking casters
- E •2-feet with electrically conductive casters
- F •1-foot with a regular swivel caster
- G •1-cheese-head screw with locking spacer
- H •1-Arm & Leg Rest Pad
- I •1 wrench

2. Slide the 5-footed plastic mid-section onto the larger end of the tripod tube. Push down on it until it comes to a stop.



Your small friend, indeed®
Health Care Logistics®
INC. 



3. Place the clamp disk with the grooves facing downward onto the 5-footed mid-section.



4. Place the locking spacer into the center opening of the clamp disk. The teeth should fit against the disk.

© Health Care Logistics, Inc. 2001

HCL-INTL is a registered trademark of Health Care Logistics, Inc.

Sheet #10656
REV.062801


Call Free 1 800 848 1633 or 888 HCL INTL® (425 4685) • Fax Free 1 800 447 2923

International Call 740 477 1686 • International Fax 740 477 2923 • Web Site <http://www.hcl-intl.com> • E-mail: hcl@hcl-intl.com


Instructions for the Arm and Leg Stand (#10656)




5. Thread the cheese-head screw into the locking spacer leaving it loose so the feet can still be easily inserted into the 5-footed plastic mid-section.



6. (A) Lift the clamp disk up and slide one of the feet with electrically conductive casters into the 5-footed plastic mid-section. (B) Align the groove on the clamp disk with the hole in the foot and push the clamp disk down until they snap together.



7. Attach the second foot with an electrically conductive caster as shown above, opposite the first.



8. Attach the two locking casters opposite one another. Attach the remaining foot with the regular caster onto the 5-footed mid-section.



Call Free 1 800 848 1633 or 888 HCL INTL® (425 4685) • Fax Free 1 800 447 2923
International Call 740 477 1686 • International Fax 740 477 2923 • Web Site <http://www.hcl-intl.com> • E-mail: hcl@hcl-intl.com

Instructions for the Arm and Leg Stand (#10656)



9. The grooves of the clamp disk must fit in the holes of the feet so that each foot is locked into the 5-footed plastic mid-section before tightening.



10. After all five feet have been placed into the plastic five footed mid-section, tighten the cheese-head screw by turning it clockwise with the wrench provided.



11. Turn the Arm and Leg Stand over and set the casters on the floor. Loosen the height adjuster by turning it counterclockwise to adjust the telescoping section to the desired height. Tighten the height adjuster by turning it clockwise to secure the telescoping section at the desired height.



12. Slide the arm and leg rest pad onto the top of the telescoping section. Push it down as far as it will go.

Your small friend, indeed®
Health Care Logistics®
INC. 

Call Free 1 800 848 1633 or 888 HCL INTL® (425 4685) • Fax Free 1 800 447 2923

International Call 740 477 1686 • International Fax 740 477 2923 • Web Site <http://www.hcl-intl.com> • E-mail: hcl@hcl-intl.com

Instructions for the Arm and Leg Stand (#10656)



13. Secure the pad to the telescoping section by turning the upper knob clockwise.



14. Adjust the Arm and Leg Rest Pad by turning the Arm and Leg Rest Adjuster counter clockwise to loosen and clockwise to tighten.



15. The Arm and Leg Stand is now ready for use.

Your small friend, indeed®
Health Care Logistics®
INC. 

Call Free 1 800 848 1633 or 888 HCL INTL® (425 4685) • Fax Free 1 800 447 2923
International Call 740 477 1686 • International Fax 740 477 2923 • Web Site <http://www.hcl-intl.com> • E-mail: hcl@hcl-intl.com