

Instructions for the Easy-Grip Tablet Crusher #7336



1. Remove the crusher lid by turning it counter clockwise. Place tablet(s) into the base of the crusher.



2. Put the crusher lid back on the crusher base turning it clockwise until the lid pressure crushes the tablet(s).



3. Rotate the lid back and forth to further crush the tablet(s).



4. The medication is ready to dispense.



Call Free: 1 800 848 1633 or 1 888 HCL INTL® (425 4685) Fax Free: 1 800 447 2923
Web Site: <http://www.hcl-intl.com> • E-mail: hcl@hcl-intl.com

Instructions for the Easy-Grip Tablet Crusher #7336



1. Remove the crusher lid by turning it counter clockwise. Place tablet(s) into the base of the crusher.



2. Put the crusher lid back on the crusher base turning it clockwise until the lid pressure crushes the tablet(s).



3. Rotate the lid back and forth to further crush the tablet(s).



4. The medication is ready to dispense.



Call Free: 1 800 848 1633 or 1 888 HCL INTL® (425 4685) Fax Free: 1 800 447 2923
Web Site: <http://www.hcl-intl.com> • E-mail: hcl@hcl-intl.com