

ASSEMBLY INSTRUCTIONS FOR


Grip 'N' Crush Crusher (#18793)



Step 1. Place the blister label side down on the countertop.



Step 2. Align the Grip 'N' Crush parallel with the long side of the tablet. This prevents the tablet from popping out from under the crusher. Push straight down to crack the tablet.



Step 3. If you encounter a stubborn tablet, place one hand on top of the crusher and press straight down.



Step 4. Rock back and forth until the tablets are crushed to the desired consistency.



Step 5. To crush multiple uncoated tablets at once, place inside a 4 or 6-mil Reclosable Bag then close to prevent dispersion of the tablet powder and keep the area clean.



Step 6. Push straight down to crack the tablets.

Call Free 1.800.848.1633 • Fax Free 1.800.447.2923

Web: GoHCL.com • Email: hcl@GoHCL.com

ASSEMBLY INSTRUCTIONS FOR

Grip 'N' Crush Crusher (#18793)



Step 7. Rock the crusher over the tablets until they are pulverized to your satisfaction.



Step 8. Then pour the crushed and pulverized tablets into the liquid or food.”

Call Free 1.800.848.1633 • Fax Free 1.800.447.2923
Web: GoHCL.com • Email: hcl@GoHCL.com