

# CLEANING GUIDE FOR

## IV Lock Boxes (18638, 18634, 18855, 18875, 20476, 20478, 20479)



18638 Double Lock-to-Pole  
IV Lock Box, Key Lock



18634 Double Lock-to-Pole  
IV Lock Box, KEDL



18855 Extended Double  
Lock-to-Pole IV Lock Box,  
Key Lock



18875 Deep Extended Double  
Lock-to-Pole IV Lock Box,  
Key Lock



20476 The Pod 2.0  
IV Lock Box  
w/Key Lock, Clear



20477 The Pod 2.0  
IV Lock Box  
w/Key Lock, Amber



20478 The Pod 2.0  
IV Lock Box  
w/KEDL Lock, Clear



20479 The Pod 2.0  
IV Lock Box  
w/KEDL Lock, Amber

Wipe surfaces with the disinfectant listed below. Allow to air dry per dwell time provided on disinfectant label.

### Suggested Disinfectant:

19954 PREempt™ One-Step Surface Cleaner and Disinfectant Wipes, 6 x 7, Canister

19954-31 PREempt™ One-Step Surface Cleaner and Disinfectant Wipes, 6 x 7, Canister, Case

Call Free: 1.800.848.1633 • Fax Free: 1.800.447.2923

GoHCL.com • hcl@GoHCL.com

© Health Care Logistics®, Inc. 2024