

INSTRUCTIONS FOR

Self-Righting Tip Caps for Exactamed™ Oral Dispensers, Blue (#7888)

- eng** EXACTAMED Oral Dispenser Tip Cap
- fre** Bouchon pour distributeur-doseur à usage oral EXACTAMED
- ger** EXACTAMED Spitzenkappe für Orale Dosiergeräte
- ita** Cappuccio per Punte da usare con il Distributore-Dosatore EXACTAMED per somministrazioni orali
- spa** Tapa de punta del surtidor oral EXACTAMED
- dut** Bescherm dopje EXACTAMED orale dispenser
- swe** EXACTAMED skyddshätta för oral dispenseringspipett



EXACTAMED Oral Dispenser Tip Cap

Cap EXACTAMED Oral Liquid Dispensers.

Specifications:

- Cap: blue or amber high-density polyethylene (HDPE)
- Self-righting
- Non-sterile
- Made in the USA



Contraindications:

None when used as described in this insert. Tip caps are a choking hazard. Keep tip caps out of the reach of children.

	Do not reuse. Product reuse can be unsafe because of the risk of: cross-infection, endotoxin reaction, patient injury due to material failure and/or chemical reaction due to cleaning or sterilization agents.		See Instructions for Use – IFU and Manuals
	Discard if packaging is damaged.		No natural rubber latex components

Call Free 1.800.848.1633 • Fax Free 1.800.447.2923

Web: GoHCL.com • Email: hcl@GoHCL.com