

# **Certification of the USP Moisture Permeation Test**

### Prepared for:

Health Care Logistics 450 Town Street Circleville, Ohio 43113-0025

## **Objective:**

To quantify moisture permeation, according to the USP Standards for the Green Deep Round Blister HCL #7025 - Lot #: 050707 and Labels for Round Blisters HCL #6125 - Lot #: 703801 of Health Care Logistics.

## **Summary:**

Ten sample blisters were packaged with dried desiccant tablets and sealed using moderate force with a sealing roller on a foam tray with pre-formed compartments. Ten empty blisters served as the control. All blisters were placed in a temperature and humidity controlled environment; taken out to equalize to room temperature and weighed on an electronic balance per USP Testing Standard Method I from Section <671> Containers-Permeation. Testing was performed for a 30-day duration. The following formula for moisture permeation was used to calculate the mg gained per day:

Moisture Permeation =  $(1/N) [W_f - W_i) - (C_f - C_i)]$ 

#### Results:

For the 30 day testing program, this package configuration was effective in preventing moisture permeation and meets the parameters provided by the USP to warrant a "Class A" package classification.

Tested by:

DevelopMed Ltd.

Kenneth F. Bober, R.Ph.

Kenneth F. Bole -

President